

# SCRUB PROGRAM CHECKLIST

What to look for in an effective scrub program.



## CERTAIN COMPONENTS OF AORN'S GUIDELINE FOR SURGICAL ATTIRE

## HOW ARAMARK CAN HELP YOU COMPLY

### Recommendation I

- Contaminated scrub attire must be bagged or containerized at facility where it was used
- Prevent contamination of laundered surgical attire during transport to the health care facility
- Surgical attire is transported in a clean vehicle and enclosed carts or containers
- Clean surgical attire is stored in a clean, enclosed cart or cabinet
- Dispensing machine may be used

- ✓ All soiled scrubs are deposited in a designated enclosed soil container on-site
- ✓ Laundered surgical attire is stored in covered carts that are poly-lined or on covered shelves to protect the clean product
- ✓ Aramark personnel are trained to handle and transport items according to OSHA's BBP standards
- ✓ All delivery vehicles are regularly cleaned and medical apparel is packaged in a poly lined cart for transport
- ✓ Aramark offers an automated enclosed dispensing system that dispenses one scrub at a time to authorized users and compliance monitoring with real-time, web-based reports

### Recommendation II

- Scrub attire is made of low-linting material

- ✓ We offer a blended fabric (65% polyester/35% cotton) or 100% polyester - both made of low linting material
- ✓ Premium quality woven fabric
- ✓ Superior fabric and construction

### Recommendation X

- All individuals entering semi restricted and restricted areas wear clean surgical attire

- ✓ scrubBank™ dispenses clean scrubs in 10 seconds to authorized users
- ✓ Scrubs are laundered in accordance with industry standards
- ✓ Wash process features EPA-registered sanitizer

**NOTE:** This checklist contains a summary of certain relevant components of AORN's Guideline for Surgical Attire. It is not intended to be a complete checklist to cover all aspects of AORN's Guidelines and should not be used as a substitute for a complete review of the recommendations and associated internal compliance.

**800-ARAMARK**  
[aramarkuniform.com/healthcare](http://aramarkuniform.com/healthcare)